

SOP 47 | Advanced administration

Perform an authorized guided movement test

Follow the existing measured walking/stopping test sequence and inspect local results.

For Mobile · Administrator / Supervisor

Before you start An authorized setup session, supported build and safe measured test area are required. Use only an approved field study.

Open Property setup > Movement tests or Saved tests

Step 1 Open Movement tests

In property setup select Movement tests. Prepare the specified straight or square course and review the recording mode.

Step 2 Start and follow the prompts

Pass start checks and start test recording. Follow each on-screen trial's walking, turning or stopping instructions and enter actual observations when requested.

Step 3 Finish the recording

Complete the prescribed trials and finish the authorized setup recording. Check pending uploads and save status.

Step 4 Inspect Saved tests

Open Saved tests on the same phone, select the recording and review the preserved trials and comparisons.


Check your result The local saved test retains its original trials and comparison results.

If something goes wrong These tests do not prove exact indoor location and do not replace operational setup or route publication.

See the picture reference on the next page.

SOP 47 Picture reference

Perform an authorized guided movement test



Step 2 · The local saved test retains its original trials and comparison results.

Workflow illustration. Follow the complete numbered instructions on the previous page.